

## Clean-up Columns



**A good sample clean-up is another important criterion for accurate and reliable results in mycotoxin analysis. Romer offers a selection of clean-up columns in different formats:**

### **MycoSep®/MultiSep®**

- 1-step, rapid clean up
- Available in 2 formats
- Stability – 18 months
- Recoveries >90%
- Part of AOAC and CEN methods
- Unique Principle – Interferences are retained, analyte passes through the column

### **StarLine™ Immuno-Affinity Column (IAC)**

- 3-step clean up
- For a wide variety of commodities and foodstuffs
- Compliant with EC commission directive 2002/26/EC
- Stability – 18 months
- Recoveries >95%
- Repeatability <5%
- Different extraction solutions can be used