



Dia Pro[®]



ACROSS RBC
POOL CELL & EXTENDED



Across Gel® RBC

Across Gel® Pool Cell

- Across Pool Cell, is used for screening of antibodies in Anti-Human Globulin (AHG) that may be clinically important in the donor serum or plasma.
- Across Pool Cell advantages;
 - Across Pool Cell has clinically important antigen panel
 - Rh, Kell, Duffy, Kidd, Lewis, P, MNS, Luth,
 - In antibody screening tests, using AHG microtube together with Across Pool Cell provides cost advantage.

Ref. No: 810113 Across Pool Cell 1x10 ml



Across Gel® Extended Panel

- Across Identicell Extended-Across Identicell Extended P products are produced additional to Across Identi Cell-Across Identi Cell P products and helps identifying multi-antibodies available free in serum or plasma samples.
- Across Extended advantages;
 - Across Pool Cell has clinically important antigen panel
 - Rh, Kell, Duffy, Kidd, Lewis, P, MNS,Luth,
 - In addition to Across Identi Cell and Across Identi Cell P can be provided/supplied in the intended process.

Ref. No: 810114 Across Identicell Extended 4x5 ml

Ref. No: 810115 Across Identicell Extended P 4x5 ml



Merkez / HeadOffice

Barbaros Mah. Tophaneliođlu Cad. Iřık Apt No: 70/1 PK 34662 Üsküdar - İstanbul / Türkiye
Tel: +90 (216) 326 70 42 • Faks: +90 (216) 340 16 89 • e-mail: diapro@diapro.com.tr

Fabrika /Factory

İnönü Mah. Gebze Plastikçiler OSB Atatürk Blv. No:25/1 PK 41400 Gebze - Kocaeli / Türkiye
Tel: +90 (262) 751 53 34 • Faks: +90 (262) 751 53 37 • e-mail: diapro@diapro.com.tr



www.diapro.com.tr