

# across

AQUA

STERILE INHALATION WATER



DIN

Fully Designed Under German Standardization

www.diapro.com.tr


DiaPro



## STERILE INHALATION WATER

Infections in patients whose treatment continues due to chronic disease affect the quality of life of the person and increase mortality.

While sterile inhalation water provides humidification of the air, it also prevents damage to the mucosal membranes. Especially in patients using nasal cannulas with asthma, COPD, tracheostomy and intubated in the intensive care unit needs oxygen support. Therefore, weekly sterilization of humidifiers and sterility of the water used are important in reducing healthcare-associated infections, morbidity and mortality by preserving mucosal integrity.

**Inhalation solutions are for single use only.**



Product Code	Description	Content	Sterilization Method	Packaging
635032	Oxygen Humidification Bottle (Volume 350 ml. Contains Sterile Highly Purified Water) with Sterile Adapter	HPW (Highly Purified Water)	 (121°C Terminal Sterilization)	6 pieces/box
635052	Oxygen Humidification Bottle (Volume 500 ml. Contains Sterile Highly Purified Water) with Sterile Adapter	HPW (Highly Purified Water)	 (121°C Terminal Sterilization)	6 pieces/box

**Head Office** Kar Plaza İçerenköy Mah. Karaman Çiftlik Yolu  
Cad. No:47 18.Kat PK 34572 Ataşehir İstanbul / TÜRKİYE  
Phone: +90 (216) 326 70 42  
Fax : +90 (216) 340 16 89  
e-mail: diapro@diapro.com.tr

**Factory** İnönü Mah. Gebze Plastikçiler OSB Atatürk Blv. No:25/1  
PK 41400 Gebze - Kocaeli / TÜRKİYE  
Phone: +90 (262) 751 53 34  
Fax : +90 (262) 751 53 37